BacT/ALERT® MB

43-02494

DENOTES REVISED SECTION

Store BacT/ALERT® MB culture bottles and MB/BacT® Enrichment Fluid between 15-30°C. Protect from direct sunlight.

See symbol glossary at end of insert.

INTENDED USE

BacT/ALERT MB culture bottles with the addition of MB/BacT Enrichment Fluid, when used with the MB/BacT® Mycobacteria Detection System (non-shaking) and the BacT/ALERT® Microbial Detection System (shaking), is a non-selective culture medium for the qualitative culture and recovery of mycobacteria from blood specimens.

SUMMARY AND EXPLANATION

BacT/ALERT MB culture bottles in combination with the MB/BacT Enrichment Fluid are designed for the cultivation of *Mycobacterium* sp. commonly isolated from blood. The medium will support growth of other aerobic organisms, including yeast, fungi, and bacteria. This complete system includes a lytic agent (Saponin), Sodium polyanethanesulfonate (SPS) and other media supplements which eliminate the processing step, prevent clotting of blood, and enhance the growth of mycobacteria. A 3 ml-5 ml volume of blood can be inoculated directly into the culture bottle. Inoculated bottles are placed into the instrument where they are incubated (35°C-37°C) and continuously monitored for microbial growth.

PRINCIPLE OF THE TEST

The MB/BacT Mycobacteria Detection System and the BacT/ALERT Microbial Detection Systems utilize a colorimetric sensor and reflected light to monitor the presence and production of carbon dioxide (CO₂) dissolved in the culture medium.¹ If microorganisms are present in the test sample, carbon dioxide is produced as the organisms metabolize the substrates in the culture medium. When growth of the microorganisms produces CO₂, the color of the gas-permeable sensor installed in the bottom of each culture bottle changes from blue-green to yellow. The lighter color results in an increase in reflectance units as monitored by the system. Bottle reflectance is monitored and recorded by the instrument every 10 minutes.

REAGENTS

For *in vitro* diagnostic use.

WARNING: To reduce the risks of accidental exposure to infectious agents, additional precautions should be taken when using the BacT/ALERT MB culture bottle in the BacT/ALERT Microbial Detection System. It is strongly recommended that the instrument be placed in a laboratory used for the routine culture of *M. tuberculosis*. For activities involving the propagation and manipulation of *M. tuberculosis* or *Mycobacterium* sp. grown in culture, Biosafety Level 3 Practices, containment equipment, and facilities are recommended.²

CAUTION: Handle specimens and inoculated culture bottles as though capable of transmitting infectious agents. All inoculated culture bottles, specimen collection needles, and blood-drawing devices should be decontaminated according to your institution's procedures.³

BacT/ALERT® MB

BacT/ALERT MB (color-coded black) – BacT/ALERT MB sterile, disposable culture bottles contain 29 ml of media and an internal sensor that detects carbon dioxide as an indicator of microbial growth. The media formulation consists of Middlebrook 7H9 Broth (0.47% w/v), Pancreatic Digest of Casein (0.1% w/v), Glycerol (1.0% w/v), Sodium polyanethanesulfonate (0.025% w/v), in purified water. Bottles contain an atmosphere of CO₂ in oxygen and nitrogen under vacuum. The composition of the media may be adjusted to meet specific performance requirements.

CAUTION: BacT/ALERT MB culture bottles are incubated at 37°C which potentially precludes the recovery of mycobacteria requiring other incubation temperatures (e.g. *M. marinum*, *M. ulcerans*, *M. haemophilum*). The following isolates were detected as positive in the MB/BacT instrument (non-shaking) and in the BacT/ALERT instrument (shaking), with the BacT/ALERT MB culture bottle during internal studies; *M. tuberculosis*, *M. avium*, *M. intracellulare*, *M. chelonae*, *M. fortuitum*, *M. kansasii*, *M. xenopi*, *M. bovis*, *M. goodii*, and *M. simiae*.

CAUTION: It is possible that certain rare, fastidious mycobacteria will not grow or may grow slowly in the BacT/ALERT MB culture bottle growth medium. If rare, fastidious mycobacteria requiring specialized media and culture conditions are suspected, alternative methods or extended incubation time, should be considered for recovery.

CAUTION: On rare occasions organisms may be encountered that grow in the BacT/ALERT MB culture bottle growth medium but do not produce sufficient carbon dioxide to be determined positive. For example, samples from patients undergoing treatment with antimycobacterial agents may show impaired growth.

MB/BacT® Enrichment Fluid

MB/BacT Enrichment Fluid vials contain a total fill volume of 5.5 ml each. MB/BacT Enrichment Fluid consists of Bovine Serum Albumin (14.5% w/v), Sodium Chloride (2.5% w/v), Oleic Acid (0.174% w/v), Saponin (4.4% w/v), in purified water.

Additional materials required

MB/BacT Mycobacteria Detection System or BacT/ALERT Microbial Detection System
Middlebrook 7H11 or other mycobacterial agar or egg-base medium

CO₂ incubator, 37°C ± 2°C

Sterile syringes and needles

Mycobactericidal disinfectant

Blood drawing apparatus

Biological safety cabinet

Sterile, disposable gloves

Disposable gowns

Disposable masks

Appropriate biohazard waste containers for materials potentially contaminated with infectious agents

Quality Control organisms

Materials available from bioMérieux

MB/BacT Mycobacteria Detection System and BacT/ALERT Microbial Detection Systems

Storage instructions

Store BacT/ALERT MB culture bottles and MB/BacT Enrichment Fluid at room temperature (15-30°C). Protect from direct sunlight. An expiration date is printed on each bottle label. Do not use the BacT/ALERT MB culture bottles or the Enrichment Fluid vials beyond the last day of the month indicated. Bottles must be at room temperature before use.

INSTRUMENTS

Review the appropriate BacT/ALERT Microbial Detection System User Manual before use.

SPECIMEN COLLECTION AND PREPARATION

Correct specimen collection is extremely important when obtaining blood culture specimens. Refer to Manual of Clinical Microbiology³ for proper specimen collection and transport procedures. Proper skin disinfection is an essential requirement to reduce the incidence of contamination.

bioMérieux recommends that inoculated culture bottles be placed into the BacT/ALERT Microbial Detection System as soon as possible after collection. Inoculated culture bottles delayed in entry should be maintained at room temperature until they can be loaded into the incubator.

BacT/ALERT MB CULTURE BOTTLE TEST PROCEDURE

Preliminary comments and precautions

- POTENTIAL INFECTIOUS TEST SPECIMEN. Observe Universal Precautions in handling and disposing of infectious materials. Handle specimens and inoculated cultures as though capable of transmitting infectious agents.
- To minimize the potential of leakage during inoculation of specimen into culture bottles, use syringes with fixed needles.
- All inoculated culture bottles, specimen collection needles, and blood-drawing devices should be decontaminated according to your institution's procedures.³
- Always carry inoculated culture bottles in a rack or basket. In case of a spill or breakage, an aerosol of mycobacteria (including *M. tuberculosis*) or other bacteria may be produced. Follow the safety guidelines for your institution and the CDC/NIH.² With the emergence of human immunodeficiency virus/acquired immunodeficiency syndrome, mycobacteremia caused by *M. tuberculosis* has occurred more frequently and blood is now considered a potential source of transmission in the laboratory.²

- Manipulation of *M. tuberculosis* complex organisms grown in culture requires Biosafety Level 3 practices, containment equipment, and facilities.²
- Take care to prevent contamination during bottle preparation, collection of the patient sample, and inoculation of the patient sample into the BacT/ALERT MB culture bottle.
- MB/BacT Enrichment Fluid contains components which are necessary to ensure optimal growth of mycobacteria and contains a lytic agent for blood cell lysis. This product must be added to all BacT/ALERT MB culture bottles.
- MB/BacT Enrichment Fluid may be added up to 24 hours before or after inoculation of the bottle with the blood specimen. If the Enrichment Fluid is added after inoculation of blood, it must be added using appropriate Biosafety Level 3 practices.²
- These culture bottles should be utilized only by trained laboratory personnel.

Bottle preparation

NOTE: Take care to prevent contamination during both bottle preparation and inoculation of the patient sample.

Prior to use, the BacT/ALERT MB culture bottles should be examined for evidence of damage or deterioration (discoloration). Bottles exhibiting evidence of damage, leakage, or deterioration should be discarded. The media in undisturbed bottles should be clear. Do not use a bottle which contains media exhibiting turbidity, a yellow sensor, or excess gas pressure; these are signs of possible contamination. Also inspect all MB/BacT Enrichment Fluid vials for evidence of damage or contamination. Do not use MB/BacT Enrichment Fluid vials exhibiting turbidity.

- Label the culture bottle with patient information. The icons on the bottle label (⊙, #, ⊕) can be defined by the user.
- Remove plastic flip-top from culture bottle and disinfect with an alcohol swab or equivalent. Allow to air dry.
- Enrichment Fluid is provided separately and must be added to the blood culture bottle for growth of mycobacteria. This addition can be done within 24 hours of specimen inoculation, or if specimen inoculation is carried out at bedside, Enrichment Fluid can be added after receiving the inoculated vials in the laboratory.
 - Remove the plastic flip-top from the Enrichment Fluid vial and disinfect with an alcohol swab or equivalent. Allow to air dry.
 - Aseptically add 1 ml of Enrichment Fluid to each BacT/ALERT MB culture bottle.
 - Employ careful, aseptic technique to avoid contamination of the Enrichment Fluid and medium in the culture bottle.

NOTE: If the Enrichment Fluid is added following inoculation of the blood specimen, then this step must be performed in a biological safety cabinet while wearing appropriate protective clothing to comply with safety standards set forth by CDC/NIH for Biosafety Level 3 guidelines.² To avoid cross contamination, use a new syringe for each bottle containing blood.

Direct draw inoculation procedure

- Prior to inoculation, disinfect the BacT/ALERT MB culture bottle top with an alcohol swab or equivalent. Allow to air dry.
- Collect the blood using a butterfly set and the BacT/ALERT Blood Collection Adapter Cap and inoculate directly into the BacT/ALERT MB culture bottle at the patient's bedside (3-5 ml per bottle). To prevent overinoculation, monitor the blood volume intake into the culture bottle, using the 5 ml incremental markings on the bottle label.
- After inoculation, swab bottle tops with gauze soaked in 2% Amphil® or other mycobactericidal agent and allow to air dry.
- Transfer the inoculated culture bottle promptly to the testing laboratory.

Blood collection tube inoculation procedure

- Blood may also be collected in a sterile SPS or heparinized tube and inoculated into the BacT/ALERT MB culture bottle in the laboratory. No prior processing of the specimen is required. Blood collected in EDTA is unacceptable since EDTA inhibits mycobacterial growth even in trace amounts.
- Disinfect the blood collection tube top and the BacT/ALERT MB culture bottle top with an alcohol swab or equivalent. Allow the septa to air dry prior to inoculation of the blood specimen.
- Using aseptic technique, remove 3-5 mls from the collection tube with a sterile syringe and inoculate into the culture bottle.
- After inoculation, swab bottle tops with gauze soaked in 2% Amphil® or other mycobactericidal agent and allow to air dry.

INCUBATION AND MONITORING

CAUTION: General caution should be taken when subculturing positive culture bottles as they could have been overfilled or contain high gas-producing organisms. Positive culture bottles contents may be under increased internal pressure. Positive culture bottles should be transiently vented before staining or disposal to release any gas produced during microbial metabolism.

- Load inoculated BacT/ALERT MB culture bottles into the MB/BacT Mycobacteria Detection System (non-shaking) or the BacT/ALERT Microbial Detection System (shaking), following the instructions provided in the User Manual.
- After BacT/ALERT MB culture bottles are loaded into the instrument, they should remain there for at least 42 days or until designated positive.

The following procedures require the use of Biosafety Level 3 practices, containment equipment and facilities.²

- When the instrument indicates that a particular cell contains a positive bottle, remove the bottle according to procedures provided in the User's Manual.
- All bottles designated positive should be smeared and subcultured for acid-fast bacilli. In a biological safety cabinet, mix bottle contents, disinfect the bottle stopper with an alcohol swab or equivalent, and allow to air dry. Remove specimen for acid-fast staining and subculture using a syringe and needle. Disinfect bottle stopper with mycobactericidal agent after subculture. If the acid-fast smear is positive, proceed with the Mycobacteria specific identification procedures employed by your institution. If the smear is negative for acid-fast bacilli, but reveals the presence of other microorganisms, a Gram stain should be performed. If both the acid-fast smear and the Gram stain are negative, indicating a possible false positive, the bottle should be reloaded into the instrument until growth of subculture or redesignation as positive or negative. Cultures which were initially determined false positive and were redesignated positive should be smeared and subcultured.
- If non-mycobacterial organisms are seen on the Gram stain, obtain another specimen for culture, if desired. The possibility of bacterial septicemia should also be considered.
- Negative cultures may be checked by smear and/or subcultured at some point prior to discarding as negative.
- Do not reuse BacT/ALERT culture bottles.** Dispose of inoculated BacT/ALERT culture bottles according to your laboratory protocol. Autoclaving and/or incinerating inoculated BacT/ALERT bottles is appropriate.²

REPORTING RESULTS

Results	Report
Signal +, AFB Smear +	AFB positive; identification pending
Signal +, AFB Smear -, Gram Stain +	No report; or specimen contaminated with non-mycobacterial organisms; unable to determine presence/absence of AFB
Signal +, AFB Smear -, Gram Stain -	No report. Subculture to determine possible false positive

QUALITY CONTROL

A Certificate of Conformance is provided with each case of BacT/ALERT MB culture bottles indicating satisfactory growth performance of *M. avium* ATCC 25291 and *M. intracellulare* ATCC 13950. Upon receipt, new lots or shipments of BacT/ALERT MB culture bottles or reagents may be tested for quality control. Refer to CLSI/NCCLS M22-A3 for appropriate Quality Control organisms. Follow the Bottle preparation and inoculation procedure:

- Add 1.0 ml MB/BacT Enrichment Fluid to each BacT/ALERT MB culture bottle required for testing.
- Inoculate representative bottles with 5.0 ml of the control organisms, diluted to 10³ CFU/ml in sterile, defibrinated sheep blood. Subculture on Middlebrook 7H11 plates to determine viability and count.
- After obtaining anticipated results, use remaining bottles for testing clinical specimens. If expected results are not obtained, contact bioMérieux Customer Service.

LIMITATIONS OF THE TEST

Many variables involved in mycobacterial testing cannot be practically controlled to provide total confidence that results obtained are due solely to proper or improper performance of any culture medium or detection system.

- Cultures obtained for primary diagnosis after initiation of therapy may produce negative results.
- Recovery of mycobacteria in the BacT/ALERT MB culture bottles is dependent on the quality of the specimen collected and the numbers of culturable organisms in the specimen volume. Since media is non-selective and does not contain antibiotics, any microorganisms present may grow in the media.
- Infrequently, BacT/ALERT MB culture bottle positives may occur due to a very high number of white blood cells being present in the blood specimen. This will result in smear and subculture negative samples.
- Organisms are often few in number. Collection of specimens over three consecutive days is generally recommended.

- Colony morphology and pigmentation should only be determined on solid media.
- BacT/ALERT MB culture bottles with a positive signal may contain one or more species of mycobacteria and/or other non-mycobacterial species. In mixed mycobacterial cultures in broth, a slow grower may be out-competed by a fast grower, or the culture may have been discarded early after the rapid grower is detected. Identification of mycobacteria present requires subculture and additional procedures to identify organisms present. The consistency of microscopic morphology in the BacT/ALERT culture bottle has not been established.
- Mycobacteria may vary in acid-fastness depending on strain, age of culture and other variables. All bottles with a positive signal from the instrument should be subcultured to both mycobacterial selective and non-selective media. Non-mycobacteria may overgrow mycobacteria present.

PERFORMANCE CHARACTERISTICS – MB/BacT (NON-SHAKING)

In-house seeded studies were performed in an MB/BacT (non-shaking) instrument using the following test organisms at levels of ≤ 10 CFU/bottle and ≤ 100 CFU/bottle in human blood from a healthy adult population.

Results of BacT/ALERT MB in MB/BacT Instrument (non-shaking)

Mycobacterium	Inoculum (CFU/bottle)	Time to Detection (days) ¹ BacT/ALERT® MB
TB Complex		
<i>M. tuberculosis</i> (4 strains) CDC 18293, CDC 18280	≤ 100	13.3 - 27.2
CDC 18283, CDC 18286	≤ 10	14.3 - 30.0 ²
<i>M. bovis</i> (2 strains) CDC 423, CDC 8131	≤ 100	14.9 - 23.8
	≤ 10	16.2 - 25.5
Photochromogens (Runyoun I)		
<i>M. kansasii</i> CDC 834, <i>M. simiae</i> CDC 524	≤ 100	10.0 - 14.6
	≤ 10	13.1 - 39.9
Scotochromogens (Runyoun II)		
<i>M. goodii</i> ATCC 14470, <i>M. xenopi</i> CDC 725	≤ 100	14.4 - 30.5
	≤ 10	17.2 - 32.6 ²
Nonchromogens (Runyoun III)		
<i>M. avium</i> (2 strains) CDC 569, CDC 1155	≤ 100	10.4 - 12.5
	≤ 10	13.8 - 14.8
<i>M. avium</i> (2 strains) CDC 440, CDC 768	≤ 100	9.2 - 13.6 ³
	≤ 10	10.9 ³ - 15.7
<i>M. intracellulare</i> (4 strains) CDC 403, VA 5, CDC 644, CDC 1406	≤ 100	8.5 - 12.1
	≤ 10	11.0 - 12.9
Rapid Growers (Runyoun IV)		
<i>M. chelonae</i> , <i>M. fortuitum</i>	≤ 100	2.2 - 2.3
	≤ 10	2.6 - 3.2

¹ Each organism was tested in triplicate and averages obtained except when noted. Values given are averages of the replicates.

² One of three bottles showed no growth; value given is average of two bottles.

³ One of three bottles contaminated; value given is average of two bottles.

PERFORMANCE CHARACTERISTICS – BacT/ALERT (SHAKING)

In-house seeded studies were performed using the following test organisms at levels of ≤ 10 CFU/bottle and ≤ 100 CFU/bottle in human blood from a healthy adult population.

Results of BacT/ALERT MB in BacT/ALERT Instrument (Shaking)

Mycobacterium	Inoculum (CFU/bottle)	Time to Detection (days) ¹ BacT/ALERT® MB
TB Complex		
<i>M. tuberculosis</i> (3 strains) CDC 18292, CDC 18280, CDC 18286	≤ 100	14.1 - 16.7
	≤ 10	14.6 - 18.4
<i>M. tuberculosis</i> CDC 18283	≤ 100	10.6
	≤ 10	13.2 ²
<i>M. bovis</i> CDC 423	≤ 100	19.8
	≤ 10	23.2
<i>M. bovis</i> CDC 8131	≤ 100	19.7
	≤ 10	NG ⁴
Photochromogens (Runyoun I)		
<i>M. kansasii</i> CDC 834, <i>M. simiae</i> CDC 524	≤ 100	7.9 - 14.4 ²
	≤ 10	8.9 - 16.8 ²
Scotochromogens (Runyoun II)		
<i>M. goodii</i> ATCC 14470, <i>M. xenopi</i> CDC 725	≤ 100	17.8 - 37.2
	≤ 10	21.3 - 43.2 ³
Nonchromogens (Runyoun III)		
<i>M. avium</i> (3 strains) CDC 569, CDC 768, CDC 1155	≤ 100	7.5 - 11.5
	≤ 10	8.5 - 11.8
<i>M. avium</i> CDC 440	≤ 100	8.4
	≤ 10	9.3 ²
<i>M. intracellulare</i> (4 strains) CDC 403, VA 5, CDC 644, CDC 1406	≤ 100	9.0 - 11.6
	≤ 10	9.9 - 13.1
Rapid Growers (Runyoun IV)		
<i>M. chelonae</i> , <i>M. fortuitum</i>	≤ 100	2.1
	≤ 10	2.3

¹ Each organism was tested in triplicate and averages obtained. Values given are a range of these averages, except *M. bovis*, *M. tuberculosis* CDC 18283, and *M. avium* CDC 440, where only the average is given.

² One of three bottles contaminated; value given is average of two bottles.

³ One of three bottles showed no growth; value given is average of two bottles.

NOTE: A list of rare and fastidious organisms recovered with BacT/ALERT culture bottles is available upon request from bioMérieux.⁴

REFERENCES

- Thorpe TC, Wilson ML, Turner JE, et al: BacT/ALERT: an automated colorimetric microbial detection system. *J Clin Micro* 1990; 28(7), 1608-1612.
- Biosafety in Microbiological and Biomedical Laboratories (BMBL) 5th Edition. U.S. Department of Health and Human Services Centers for Disease Control and Prevention and National Institutes of Health. Fifth Edition. US Government Printing Office. Washington: Feb 2007.
- Metchock B, Nolte FS, Wallace Jr RJ: Mycobacterium, in Murray PR, (ed): *Manual of Clinical Microbiology*, ed 7. Washington, DC, American Society for Microbiology, 1999, pp 399-437.
- Rare Organism Club, bioMérieux, Inc., 2007.

AVAILABILITY

bioMérieux

BacT/ALERT® MB 25/case  251011

MB/BacT® Enrichment Fluid 5 vials/kit  259877

For technical assistance in the U.S.A., contact bioMérieux Customer Service at 1-800-682-2666. Outside the U.S.A., contact your local bioMérieux Representative.

bioMérieux, the blue logo, BacT/ALERT, and MB/BacT are used, pending, and/or registered trademarks belonging to bioMérieux SA or one of its subsidiaries.

Amphyl is a registered trademark of Linden Corporation.

ATCC is a registered trademark of the American Type Culture Collection.

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1. Thorpe TC, Wilson ML, Turner JE, et al: BacT/ALERT: an automated colorimetric microbial detection system. *J Clin Micro* 1990; 28(7), 1608-1612.
2. Biosafety in Microbiological and Biomedical Laboratories (BMBL) 5th Edition. U.S. Department of Health and Human Services Centers for Disease Control and Prevention and National Institutes of Health. Fifth Edition. US Government Printing Office. Washington: Feb 2007.
3. Metchock B, Nolte FS, Wallace Jr RJ: Mycobacterium, in Murray PR, (ed): *Manual of Clinical Microbiology*, ed 7. Washington, DC, American Society for Microbiology, 1999, pp 399-437.
4. Rare Organism Club, bioMérieux, Inc., 2007.

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- Thorpe TC, Wilson ML, Turner JE, et al: BacT/ALERT: an automated colorimetric microbial detection system. *J Clin Micro* 1990; 28(7): 1608-1612.
- Biosafety in Microbiological and Biomedical Laboratories (BMBL) 5th Edition. U.S. Department of Health and Human Services Centers for Disease Control and Prevention and National Institutes of Health. Fifth Edition. US Government Printing Office. Washington: Feb 2007.
- Metchock B, Nolte FS, Wallace Jr RJ: Mycobacterium, in Murray PR, (ed): *Manual of Clinical Microbiology*, ed 7. Washington, DC, American Society for Microbiology, 1999, pp 399-437.
- Rare Organism Club, bioMérieux, Inc., 2007.

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REFERANSELISTE

1. Thorpe TC, Wilson ML, Turner JE, et al: BacT/ALERT: an automated colorimetric microbial detection system. *J Clin Micro* 1990; 28(7), 1608-1612.
2. Biosafety in Microbiological and Biomedical Laboratories (BMBL) 5th Edition. U.S. Department of Health and Human Services Centers for Disease Control and Prevention and National Institutes of Health. Fifth Edition. US Government Printing Office. Washington: Feb 2007.
3. Metchnik B, Nolte FS, Wallace Jr RJ: Mycobacterium, in Murray PR, (ed): *Manual of Clinical Microbiology*, ed 7. Washington, DC, American Society for Microbiology, 1999, pp 399-437.
4. Rare Organism Club, bioMérieux, Inc., 2007.

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1. Thorpe TC, Wilson ML, Turner JE, et al: BacT/ALERT: an automated colorimetric microbial detection system. *J Clin Micro* 1990; 28(7), 1608-1612.
2. Biosafety in Microbiological and Biomedical Laboratories (BMBL) 5th Edition. U.S. Department of Health and Human Services Centers for Disease Control and Prevention and National Institutes of Health. Fifth Edition. US Government Printing Office. Washington: Feb 2007.
3. Metchock B, Nolte FS, Wallace Jr RJ: Mycobacterium, in Murray PR, (ed): *Manual of Clinical Microbiology*, ed 7. Washington, DC: American Society for Microbiology, 1999, pp 399-437.
4. Rare Organism Club, bioMérieux, Inc., 2007.

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©BIOMÉRIEUX 2000, 2004, 2007, 2008
Amphyl är ett registrerat varumärke som tillhör Linden Corporation.
bioMérieux, S.A.