## Instrument

A Bact/ALERT Reflectance Standard kit is provided with each instrument for the QC and Calibration procedures. All quality control should be part of normal system maintenance. Refer to the Bact/ALERT User Manual for more information.

## LIMITATIONS OF THE TEST

Many variables involved in mycobacterial testing cannot be practically controlled to provide total confidence that results obtained are due solely to proper or improper performance of any culture medium or detection system.

- Cultures obtained for primary diagnosis after initiation of antimicrobial therapy may produce negative results.
- Recovery of mycobacteria in the Bact/ALERT MP bottles is dependent on the quality of the specimen collected, the numbers of culturable organisms in the specimen volume, and the method of processing. Adherence to procedural instructions is critical for optimum recovery of mycobacteria.
- Organisms are often few in numbers. Collection of specimens over three consecutive days is generally recommended.
- Improper processing/decontamination or delays in transport of specimens can result in an overgrowth of bacteria which could lead to a contaminated broth culture. The entire culture medium can become contaminated if even one viable bacterium is present in the specimen after processing.
- Contaminants resistant to the antibiotics contained in the MB/Bact Antibiotic Supplement may be encountered.
- Colony morphology and pigmentation should only be determined on solid media.
- Bact/ALERT MP culture bottles with a positive signal may contain one or more species of mycobacteria and/or other non-mycobacterial species. In mixed mycobacteria cultures in broth, a slow grower may be out-competed by a fast grower or the culture discarded early after the rapid grower is detected. Identification of mycobacteria present requires subculture and additional procedures to identify organisms present. The consistency of microscopic morphology in Bact/ALERT MP culture bottles has not been established.
- Contamination with saprophytic mycobacteria in tap water or other laboratory reagents and equipment may cause "clinically" false positive results (recovery of mycobacteria not in the clinical specimen).
- Decontamination by the N-Acetyl-L-Cysteine-Sodium Hydroxide method is recommended. Other decontamination methods have not been tested in conjunction with Bact/ALERT MP culture medium. Digestant-decontaminant reagents may have harmful effects on mycobacteria.
- Mycobacteria may vary in acid-fastness depending on strain, age of culture, and other variables. All bottles with a positive Bact/ALERT MP signal or appearing turbid should be subcultured to both mycobacterial selective and nonselective media. Non-mycobacterial species may overgrow mycobacteria present. Such culture bottles should be decontaminated and recultured.
- Bact/ALERT MP culture bottles are incubated at 35°C precluding the recovery of mycobacteria requiring other incubation temperatures (e.g., *M. marinum*, *M. ulcerans*, and *M. haemophilum*). Recovery of such organisms requires additional culture methods. Organisms with special growth requirements (e.g., *M. haemophilum*) may not be recovered in Bact/ALERT MP bottles when incubated at the appropriate temperature. The following isolates have been recovered in analytical studies: *M. tuberculosis*, *M. bovis*, *M. avium-intracellulare*, *M. malmoense*, *M. kansasii*, *M. simiae*, *M. fortuitum*, *M. chelonae*, *M. gordonae*, and *M. xenopi*. The ability of the Bact/ALERT MP bottles to recover other uncommon isolates is unknown.
- Ground tissue samples are difficult to inoculate into Bact/ALERT MP culture bottles with safety needles and syringes. In order to ensure adequate inoculation, tissue specimens must be ground finely enough to permit adequate inoculation of the material into the Bact/ALERT MP bottles. In addition, it may be advantageous to use larger gauge safety needles, or remove the removable closure, aseptically add the ground specimen and apply a Reseal to the Bact/ALERT MP culture bottle. If contamination of ground tissue samples is suspected, decontaminate and process prior to adding to Bact/ALERT MP bottles. This will allow maximum recovery of mycobacteria that may be present.
- Prolonged incubation or the addition of Mycobactin J may be required when *M. malmoense* or *M. genavense* is suspected. *M. malmoense* and *M. genavense* are known to require extended incubation past the normal schedules employed by many U.S. laboratories.<sup>2</sup>
- The Antibiotic Supplement, although necessary for the testing of most specimens, may inhibit the growth of some mycobacteria.
- False negative rates from bioMérieux's clinical evaluation were 0.2 %.

## PERFORMANCE CHARACTERISTICS OF THE TEST

Seeded studies were performed using the following organisms at levels of  $\leq 10^6$  CFU/bottle and  $\leq 10^2$  CFU/bottle.

Microorganism	Inoculum (CFU/bottle)	Time to Detection (days) <sup>1</sup>	
		Bact/ALERT MP (Plastic)	MP
<b>TB Complex</b>			
<i>M. tuberculosis</i> (4 strains)			
	$\leq 10^6$	5.6 – 10.3	
	$\leq 10^2$	13.7 – 18.8	
<i>M. bovis</i> (2 strains)			
	$\leq 10^6$	7.1 – 7.8	
	$\leq 10^2$	14.4 – 17.4	
<b>Photochromogens (Runyoun I)</b>			
<i>M. kansasii</i> (2 strains)			
	$\leq 10^6$	7.4 – 8.8	
	$\leq 10^2$	13.1 – 15.0	
<i>M. simiae</i>			
	$\leq 10^6$	4.1	
	$\leq 10^2$	8.9	
<b>Scotochromogens (Runyoun II)</b>			
<i>M. gordonae</i>			
	$\leq 10^6$	14.1	
	$\leq 10^2$	23.7	
<i>M. xenopi</i> (2 strains)			
	$\leq 10^6$	10.3 – 13.8	
	$\leq 10^2$	19.9 – 25.9 <sup>2</sup>	
<i>M. scrofulaceum</i> (2 strains)			
	$\leq 10^6$	12.6 – 23.3	
	$\leq 10^2$	20.7 – NG <sup>3</sup>	
<b>Nonchromogens (Runyoun III)</b>			
<i>M. avium</i> (4 strains)			
	$\leq 10^6$	4.2 – 6.8	
	$\leq 10^2$	10.2 – 15.9	
<i>M. intracellulare</i> (4 strains)			
	$\leq 10^6$	5.6 – 7.1	
	$\leq 10^2$	10.8 – 16.2	
<i>M. malmoense</i>			
	$\leq 10^6$	15.4	
	$\leq 10^2$	21.9	
<b>Rapid Growers (Runyoun IV)</b>			
<i>M. chelonae</i>			
	$\leq 10^6$	2.2	
	$\leq 10^2$	4.1	
<i>M. fortuitum</i> (2 strains)			
	$\leq 10^6$	3.5 – 4.1	
	$\leq 10^2$	3.9 – 8.9	

<sup>1</sup> Each organism was tested with MB/Bact Antibiotic Supplement. Testing was performed for single bottles at  $10^6$  and in triplicate at  $\leq 10^2$ . Data listed represent a range of values from multiple strains, except for *M. simiae*, *M. gordonae*, *M. malmoense*, and *M. chelonae* where data for individual organism is given.

<sup>2</sup> One of three bottles showed no growth after 42 days; value given is average of two bottles.

<sup>3</sup> One organism tested showed no growth for three of three bottles after 42 days.

## REFERENCES

- Global tuberculosis control: surveillance, planning, financing. WHO report 2007. Geneva, World Health Organization (WHO/HTM/TB/2007.376).
- Metchock BG, Nolte FS, and Wallace RJ, Jr.: Mycobacterium, in Murray PR, Baron EJ, Pfaller MA, et al (eds): *Manual of Clinical Microbiology*, ed 7. Washington, DC, American Society for Microbiology, 1999, pp 399-437.
- National Center for Health Statistics. Health, United States, 2007. With Chartbook on Trends in the Health of Americans. Hyattsville, MD: 2007.
- Eisenstadt J, Hall GS, Gibson SM, et al: Mycobacterium tuberculosis and other nontuberculosis mycobacteria, in Mahon CR, Manuselis G Jr (eds): *Textbook of Diagnostic Microbiology*. Philadelphia, PA, Saunders, 1995, pp 635-676.
- Roberts GD, Koneman EW, Kim YK: Mycobacterium in Balows A, Hausler WJ Jr, Herrmann KL, et al (eds): *Manual of Clinical Microbiology*, ed 5. Washington, DC, American Society for Microbiology, 1991, pp 304-339.
- Essential components of a tuberculosis program: recommendations of the Advisory Council for the Elimination of Tuberculosis. *MMWR* 1995;44 (No. RR-11):13.
- Tenover FC, Crawford JT, Huebner RE, et al: The resurgence of tuberculosis: is your laboratory ready? *J Clin Micro* 31(4): 767-770, 1993.
- Thorpe TC, Wilson ML, Turner JE, et al: Bact/Alert: an automated colorimetric microbial detection system. *J Clin Micro* 1990;28(7):1608-1612.
- Biosafety in Microbiological and Biomedical Laboratories (BMBL) 5th Edition. U.S. Department of Health and Human Services Centers for Disease Control and Prevention and National Institutes of Health. Fifth Edition. US Government Printing Office. Washington: Feb 2007.
- Della-Latta P (ed): Mycobacteriology and Antimycobacterial Susceptibility Testing, in Isenberg HD (ed). *Clinical Microbiology Procedures Handbook*, vol 2. Washington, DC, ASM Press; 2004 sect 7.

## AVAILABILITY

### bioMérieux

#### Bact/ALERT® MP

100/case

REF 259797

#### MB/Bact® Antibiotic Supplement Kit

100 tests/kit

REF 259760

#### MB/Bact® Antibiotic Supplement

5 vials/kit

#### MB/Bact® Reconstitution Fluid

5 vials/kit

#### Bact/ALERT® Reseal

100 units/box

REF 259787

For technical assistance in the U.S.A., contact bioMérieux Customer Service at 1-800-682-2666. Outside the U.S.A., contact your local bioMérieux Representative.

bioMérieux, the blue logo, and Bact/ALERT are used, pending, and/or registered trademarks belonging to bioMérieux SA or one of its subsidiaries.

Amphyl is a registered trademark of the Linden Corporation.

ATCC is a registered trademark of American Type Culture Collection.

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## LITTERATURENHVISNINGER

1. Global tuberculosis control: surveillance, planning, financing. WHO report 2007. Geneva, World Health Organization (WHO/HTM/TB/2007.376).
2. Metchock BG, Nolte FS, and Wallace RJ, Jr.: Mycobacterium, in Murray PR, Baron EJ, Pfaller MA, et al (eds): *Manual of Clinical Microbiology*, ed 7. Washington, DC, American Society for Microbiology, 1999, pp 399-437.
3. National Center for Health Statistics. Health, United States, 2007. With Chartbook on Trends in the Health of Americans. Hyattsville, MD: 2007.
4. Eisenstadt J, Hall GS, Gibson SM, et al: Mycobacterium tuberculosis and other nontuberculosis mycobacteria, in Mahon CR, Manuelli G Jr (eds): *Textbook of Diagnostic Microbiology*. Philadelphia, PA, Saunders, 1995, pp 635-676.
5. Roberts GD, Koneman EW, Kim YK: Mycobacterium in Balows A, Hausler WJ Jr, Herrmann KL, et al (eds): *Manual of Clinical Microbiology*, ed 5. Washington, DC, American Society for Microbiology, 1991, pp 304-339.
6. Essential components of a tuberculosis program: recommendations of the Advisory Council for the Elimination of Tuberculosis. *MMWR* 1995;44 (No. RR-11):13.
7. Tenover FC, Crawford JT, Huebner RE, et al: The resurgence of tuberculosis: is your laboratory ready? *J Clin Micro* 31 (4): 767-770, 1993.
8. Thorpe TC, Wilson ML, Turner JE, et al: BacT/Alert: an automated colorimetric microbial detection system. *J Clin Micro* 1990;28(7),1608-1612.
9. Biosafety in Microbiological and Biomedical Laboratories (BMBL) 5th Edition. U.S. Department of Health and Human Services Centers for Disease Control and Prevention and National Institutes of Health. Fifth Edition. US Government Printing Office, Washington: Feb 2007.
10. Della-Latta P (ed): Mycobacteriology and Antimycobacterial Susceptibility Testing, in Isenberg HD (ed). *Clinical Microbiology Procedures Handbook*, vol 2. Washington, DC, ASM Press; 2004 sect 7.

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## LITERATURHINWEISE

- Global tuberculosis control: surveillance, planning, financing. WHO report 2007. Geneva, World Health Organization (WHO/HTM/TB/2007.376).
- Metchock BC, Nolte FS, and Wallace RJ, Jr.: Mycobacterium, in Murray PR, Baron EJ, Pfaller MA, et al (eds): *Manual of Clinical Microbiology*, ed 7. Washington, DC, American Society for Microbiology, 1999, pp 399-437.
- National Center for Health Statistics. Health, United States, 2007. With Chartbook on Trends in the Health of Americans. Hyattsville, MD: 2007.
- Eisenstadt J, Hall GS, Gibson SM, et al: Mycobacterium tuberculosis and other nontuberculous mycobacteria, in Mahon CR, Manuells G Jr (eds): *Textbook of Diagnostic Microbiology*. Philadelphia, PA, Saunders, 1995, pp 635-676.
- Roberts GD, Koneman EW, Kim YK: Mycobacterium in Balows A, Hausler WJ Jr, Herrmann KL, et al (eds): *Manual of Clinical Microbiology*, ed 5. Washington, DC, American Society for Microbiology, 1991, pp 304-339.
- Essential components of a tuberculosis program: recommendations of the Advisory Council for the Elimination of Tuberculosis. *MMWR* 1995;44 (No. RR-11):13.
- Tenover FC, Crawford JT, Huebner RE, et al: The resurgence of tuberculosis: is your laboratory ready? *J Clin Micro* 31(4): 767-770, 1993.
- Thorpe TC, Wilson ML, Turner JE, et al: BacT/Alert: an automated colorimetric microbial detection system. *J Clin Micro* 1990;28(7): 1608-1612.
- Biosafety in Microbiological and Biomedical Laboratories (BMBL) 5th Edition. U.S. Department of Health and Human Services Centers for Disease Control and Prevention and National Institutes of Health. Fifth Edition. US Government Printing Office. Washington: Feb 2007.
- Della-Latta P (ed): Mycobacteriology and Antimycobacterial Susceptibility Testing, in ISENBERG HD (ed). *Clinical Microbiology Procedures Handbook*, vol 2. Washington, DC, ASM Press; 2004 sect 7.

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- Global tuberculosis control: surveillance, planning, financing. WHO report 2007. Geneva, World Health Organization (WHO/HTM/TB/2007.376).
- Metchock BG, Nolte FS, and Wallace RJ, Jr.: Mycobacterium, in Murray PR, Baron EJ, Pfaller MA, et al (eds): *Manual of Clinical Microbiology*, ed 7. Washington, DC, American Society for Microbiology, 1999, pp 399-437.
- National Center for Health Statistics. Health, United States, 2007. With Chartbook on Trends in the Health of Americans. Hyattsville, MD: 2007.
- Eisenstadt J, Hall GS, Gibson SM, et al: Mycobacterium tuberculosis and other nontuberculosis mycobacteria, in Mahon CR, Manuvelis G Jr (eds): *Textbook of Diagnostic Microbiology*. Philadelphia, PA, Saunders, 1995, pp 635-676.
- Roberts GD, Koneman EW, Kim YK: Mycobacterium in Balows A, Hausler WJ Jr, Herrmann KL, et al (eds): *Manual of Clinical Microbiology*, ed 5. Washington, DC, American Society for Microbiology, 1991, pp 304-339.
- Essential components of a tuberculosis program: recommendations of the Advisory Council for the Elimination of Tuberculosis. *MMWR* 1995;44 (No. RR-11):13.
- Tenover FC, Crawford JT, Huebner RE, et al: The resurgence of tuberculosis: is your laboratory ready? *J Clin Micro* 31(4): 767-770, 1993.
- Thorpe TC, Wilson ML, Turner JE, et al: BacT/Alert: an automated colorimetric microbial detection system. *J Clin Micro* 1990;28(7):1608-1612.
- Biosafety in Microbiological and Biomedical Laboratories (BMBL) 5th Edition. U.S. Department of Health and Human Services Centers for Disease Control and Prevention and National Institutes of Health. Fifth Edition. US Government Printing Office. Washington: Feb 2007.
- Della-Latta P (ed): Mycobacteriology and Antimycobacterial Susceptibility Testing, in Isenberg HD (ed). *Clinical Microbiology Procedures Handbook*, vol 2. Washington, DC, ASM Press; 2004 sect 7.

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- Global tuberculosis control: surveillance, planning, financing. WHO report 2007. Geneva, World Health Organization (WHO/HTM/TB/2007.376).
- Metchock BG, Nolte FS, and Wallace RJ, Jr.: Mycobacterium, in Murray PR, Baron EJ, Pfaller MA, et al (eds): *Manual of Clinical Microbiology*, ed 7. Washington, DC, American Society for Microbiology, 1999, pp 399-437.
- National Center for Health Statistics. Health, United States, 2007. With Chartbook on Trends in the Health of Americans. Hyattsville, MD: 2007.
- Eisenstadt J, Hall GS, Gibson SM, et al: Mycobacterium tuberculosis and other nontuberculosis mycobacteria, in Mahon CR, Manuvelis G Jr (eds): *Textbook of Diagnostic Microbiology*. Philadelphia, PA, Saunders, 1995, pp 635-676.
- Roberts GD, Koneman EW, Kim YK: Mycobacterium in Balows A, Hausler WJ Jr, Herrmann KL, et al (eds): *Manual of Clinical Microbiology*, ed 5. Washington, DC, American Society for Microbiology, 1991, pp 304-339.
- Essential components of a tuberculosis program: recommendations of the Advisory Council for the Elimination of Tuberculosis. *MMWR* 1995;44 (No. RR-11):13.
- Tenover FC, Crawford JT, Huebner RE, et al: The resurgence of tuberculosis: is your laboratory ready? *J Clin Micro* 31(4): 767-770, 1993.
- Thorpe TC, Wilson ML, Turner JE, et al: BacT/Alert: an automated colorimetric microbial detection system. *J Clin Micro* 1990;28(7):1608-1612.
- Biosafety in Microbiological and Biomedical Laboratories (BMBL) 5th Edition. U.S. Department of Health and Human Services Centers for Disease Control and Prevention and National Institutes of Health. Fifth Edition. US Government Printing Office. Washington: Feb 2007.
- Della-Latta P (ed): Mycobacteriology and Antimicrobial Susceptibility Testing, in Isenberg HD (ed). *Clinical Microbiology Procedures Handbook*, vol 2. Washington, DC, ASM Press; 2004 sect 7.

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## BIBLIOGRAFIA

1. Global tuberculosis control: surveillance, planning, financing. WHO report 2007. Geneva, World Health Organization (WHO/HTM/TB/2007.376).
2. Metchnikov BG, Nolte FS, and Wallace RJ, Jr.: Mycobacterium, in Murray PR, Baron EJ, Pfaller MA, et al (eds): *Manual of Clinical Microbiology*, ed 7. Washington, DC, American Society for Microbiology, 1999, pp 399-437.
3. National Center for Health Statistics. Health, United States, 2007. With Chartbook on Trends in the Health of Americans. Hyattsville, MD: 2007.
4. Eisenstadt J, Hall GS, Gibson SM, et al: Mycobacterium tuberculosis and other nontuberculous mycobacteria, in Mahon CR, Manuvelis G Jr (eds): *Textbook of Diagnostic Microbiology*. Philadelphia, PA, Saunders, 1995, pp 635-676.
5. Roberts GD, Koneman EW, Kim YK: Mycobacterium in Balows A, Hausler WJ Jr, Herrmann KL, et al (eds): *Manual of Clinical Microbiology*, ed 5. Washington, DC, American Society for Microbiology, 1991, pp 304-339.
6. Essential components of a tuberculosis program: recommendations of the Advisory Council for the Elimination of Tuberculosis. *MMWR* 1995;44 (No. RR-11):13.
7. Tenover FC, Crawford JT, Huebner RE, et al: The resurgence of tuberculosis: is your laboratory ready? *J Clin Micro* 31(4): 767-770, 1993.
8. Thorpe TC, Wilson ML, Turner JE, et al: BacT/Alert: an automated colorimetric microbial detection system. *J Clin Micro* 1990;28(7):1608-1612.
9. Biosafety in Microbiological and Biomedical Laboratories (BMBL) 5th Edition. U.S. Department of Health and Human Services Centers for Disease Control and Prevention and National Institutes of Health. Fifth Edition. US Government Printing Office. Washington: Feb 2007.
10. Della-Latta P (ed): Mycobacteriology and Antimicrobial Susceptibility Testing, in Isenberg HD (ed). *Clinical Microbiology Procedures Handbook*, vol 2. Washington, DC, ASM Press; 2004 sect 7.

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## REFERANSELISTE

- Global tuberculosis control: surveillance, planning, financing. WHO report 2007. Geneva, World Health Organization (WHO/HTM/TB/2007.376).
- Metchock BG, Nolte FS, and Wallace RJ, Jr.: Mycobacterium, in Murray PR, Baron EJ, Pfaller MA, et al (eds): *Manual of Clinical Microbiology*, ed 7. Washington, DC, American Society for Microbiology, 1999, pp 399-437.
- National Center for Health Statistics. Health, United States, 2007. With Chartbook on Trends in the Health of Americans. Hyattsville, MD: 2007.
- Eisenstadt J, Hall GS, Gibson SM, et al: Mycobacterium tuberculosis and other nontuberculous mycobacteria, in Mahon CR, Manuvelis G Jr (eds): *Textbook of Diagnostic Microbiology*. Philadelphia, PA, Saunders, 1995, pp 635-676.
- Roberts GD, Koneman EW, Kim YK: Mycobacterium in Balows A, Hausler WJ Jr, Herrmann KL, et al (eds): *Manual of Clinical Microbiology*, ed 5. Washington, DC, American Society for Microbiology, 1991, pp 304-339.
- Essential components of a tuberculosis program: recommendations of the Advisory Council for the Elimination of Tuberculosis. *MMWR* 1995;44 (No. RR-11):13.
- Tenover FC, Crawford JT, Huebner RE, et al: The resurgence of tuberculosis: is your laboratory ready? *J Clin Micro* 31(4): 767-770, 1993.
- Thorpe IC, Wilson ML, Turner JE, et al: Bact/Alert: an automated colorimetric microbial detection system. *J Clin Micro* 1990;28(7), 1608-1612.
- Biosafety in Microbiological and Biomedical Laboratories (BMBL) 5th Edition. U.S. Department of Health and Human Services Centers for Disease Control and Prevention and National Institutes of Health. Fifth Edition. US Government Printing Office. Washington: Feb 2007.
- Della-Latta P (ed): Mycobacteriology and Antimycobacterial Susceptibility Testing, in Isenberg HD (ed). *Clinical Microbiology Procedures Handbook*, vol 2. Washington, DC, ASM Press; 2004 sect 7.

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- Global tuberculosis control: surveillance, planning, financing. WHO report 2007. Geneva, World Health Organization (WHO)/HTM/TB/2007.376).
- Metchock BG, Nolte FS, and Wallace RJ, Jr.: Mycobacterium, in Murray PR, Baron EJ, Pfaller MA, et al (eds): *Manual of Clinical Microbiology*, ed 7. Washington, DC, American Society for Microbiology, 1999, pp 399-437.
- National Center for Health Statistics. Health, United States, 2007. With Chartbook on Trends in the Health of Americans. Hyattsville, MD: 2007.
- Eisenstadt J, Hall GS, Gibson SM, et al: Mycobacterium tuberculosis and other nontuberculosis mycobacteria, in Mahon CR, Manuselis G Jr (eds): *Textbook of Diagnostic Microbiology*. Philadelphia, PA, Saunders, 1995, pp 635-676.
- Roberts GD, Koneman EW, Kim YK: Mycobacterium in Balows A, Hausler WJ Jr, Herrmann KL, et al (eds): *Manual of Clinical Microbiology*, ed 5. Washington, DC, American Society for Microbiology, 1991, pp 304-339.
- Essential components of a tuberculosis program: recommendations of the Advisory Council for the Elimination of Tuberculosis. *MMWR* 1995;44 (No. RR-11):13.
- Tenover FC, Crawford JT, Huebner RE, et al: The resurgence of tuberculosis: is your laboratory ready? *J Clin Micro* 31(4): 767-770, 1993.
- Thorpe TC, Wilson ML, Turner JE, et al: BacT/Alert: an automated colorimetric microbial detection system. *J Clin Micro* 1990;28(7), 1608-1612.
- Biosafety in Microbiological and Biomedical Laboratories (BMBL) 5th Edition. U.S. Department of Health and Human Services Centers for Disease Control and Prevention and National Institutes of Health. Fifth Edition. US Government Printing Office. Washington: Feb 2007.
- Della-Latta P (ed): Mycobacteriology and Antimycobacterial Susceptibility Testing, in Isenberg HD (ed). *Clinical Microbiology Procedures Handbook*, vol 2. Washington, DC, ASM Press; 2004 sect 7.

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
## BIBLIOGRAFIA

- Global tuberculosis control: surveillance, planning, financing. WHO report 2007. Geneva, World Health Organization (WHO/HTM/TB/2007.376).
- Metchock BG, Nolte FS, and Wallace RJ, Jr.: Mycobacterium, in Murray PR, Baron EJ, Pfaller MA, et al (eds): *Manual of Clinical Microbiology*, ed 7. Washington, DC, American Society for Microbiology, 1999, pp 399-437.
- National Center for Health Statistics. Health, United States, 2007. With Chartbook on Trends in the Health of Americans. Hyattsville, MD: 2007.
- Eisenstadt J, Hall GS, Gibson SM, et al: Mycobacterium tuberculosis and other nontuberculosis mycobacteria, in Mahon CR, Manuvelis G Jr (eds): *Textbook of Diagnostic Microbiology*. Philadelphia, PA, Saunders, 1995, pp 635-676.
- Roberts GD, Koneman EW, Kim YK: Mycobacterium in Balows A, Hausler WJ Jr, Herrmann KL, et al (eds): *Manual of Clinical Microbiology*, ed 5. Washington, DC, American Society for Microbiology, 1991, pp 304-339.
- Essential components of a tuberculosis program: recommendations of the Advisory Council for the Elimination of Tuberculosis. *MMWR* 1995;44 (No. RR-11):13.
- Tenover FC, Crawford JT, Huebner RE, et al: The resurgence of tuberculosis: is your laboratory ready? *J Clin Micro* 31(4): 767-770, 1993.
- Thorpe TC, Wilson ML, Turner JE, et al: BacT/Alert: an automated colorimetric microbial detection system. *J Clin Micro* 1990;28(7), 1608-1612.
- Biosafety in Microbiological and Biomedical Laboratories (BMBL) 5th Edition. U.S. Department of Health and Human Services Centers for Disease Control and Prevention and National Institutes of Health. Fifth Edition. US Government Printing Office. Washington: Feb 2007.
- Della-Latta P (ed): Mycobacteriology and Antimycobacterial Susceptibility Testing, in Isenberg HD (ed). *Clinical Microbiology Procedures Handbook*, vol 2. Washington, DC, ASM Press; 2004 sect 7.

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©BIOMÉRIEUX 2000, 2004, 2007, 2008
Amphyl är ett registrerat varumärke som tillhör Linden Corporation.
bioMérieux, Inc.
bioMérieux, S.A.